

Valid until 30/06/2023



OFFERS 2023



DISCOVER THE BEST OF
SWISS TECHNOLOGY

SWISS  MADE

NEW CONTRA-ANGLE

NOVA

Dear Practitioners,

Bien-Air has been accompanying you at the heart of your profession for over 60 years. As a Swiss family-owned company with an international reputation, we are committed to developing innovative, reliable and easy-to-use rotary instruments and systems year after year.

Beyond the quality of our products, we offer a complete range of services to make your customer experience a success.

In your practice or remotely, we are there for you every day.

From product demonstrations to the installation of our systems. From advice on maintenance and servicing of our instruments, to training of your assistants, our teams and dealers network will help you find THE right solution for your needs.

Benefit from a unique and innovative customer experience with our Repair Center which will answer all your technical questions and guarantee certified repairs by our Bien-Air technicians.

At Bien-Air you are not only investing in a product but also in a quality of service. We aim to provide you with a first class customer service. Whatever your question, you will always find your dedicated contact person.

The Bien-Air UK team wishes you the best for this year 2023

OUR TEAM AT YOUR SERVICE



+ 44 (0)1293 550 200
ba-uk@bienair.com



Carl Wood

Area Sales Manager South
07967 593 655
carl.wood@bienair.com

Paul Taylor

Area Sales Manager North
07553 106 005
paul.taylor@bienair.com



Scan the QR code or
visit www.bienair.com

NEW
TORNADO

NEW
PROTECTIVE SHIELD
4X HARDER

SWISS  MADE

TURBINES



Power is control



Extended service life


Tornado are now the most powerful turbines of the brand thanks to its **Steady Torque** technology, which combines optimal air pressure distribution and exhaust flow.

The **wear-resistant ceramic ball** bearings guarantee an extended service life. Every Bien-Air turbine is individually, laser-balanced to perform at very high speeds with virtually zero vibration.



Increased reliability & precision

The **Accu-Chuck PreciPlus** chuck design ensures the most secure clamping of the bur. Noise is significantly reduced, thereby lending the Tornado unequalled practitioner and patient comfort.



Anti-Return Valve

Reduced risk of cross-contamination thanks to the non-return valve and **Sealed Head** protection. These reduce the aspiration of fluids into the head of the instrument.



TORNADO
£ 725_{+VAT}

1 x TORNADO LED
REF. 1600908-001
1 x TORNADO LK
REF. 1600907-001


- Protective Shield (PVD Treatment)
- More powerful - 30 Watts
- Quieter - 55 dBA
- Head size - Ø 12,2 mm, h 13 mm
- 4 asymmetrical mixed air/water sprays
- LED lighting



TORNADO^S
£ 725_{+VAT}

1 x TORNADO^S LED
REF. 1600933-001
1 x TORNADO^S LK
REF. 1600932-001

- Protective Shield (PVD Treatment)
- More powerful - 27 Watts
- Quieter - 55 dBA
- Head size - Ø 10,8 mm, h 12,1 mm
- 4 asymmetrical mixed air/water sprays
- LED lighting



DUO PACK TORNADO PVD
+ FREE COUPLING
(CHOOSE BETWEEN UNIFIX OR MULTIFLEX)

£ 1 350_{+VAT}

2 x TORNADO LED
REF. 1600908-001
OR 2 x TORNADO LK
REF. 1600907-001

+ COUPLING LK
REF. 1600902-001
OR COUPLING UNIFIX 4 HOLE
REF. 1600363-001



DUO PACK TORNADO^S PVD
+ FREE COUPLING
(CHOOSE BETWEEN UNIFIX OR MULTIFLEX)

£ 1 350_{+VAT}

2 x TORNADO^S LED
REF. 1600933-001
OR 2 x TORNADO^S LK
REF. 1600932-001

+ COUPLING LK
REF. 1600902-001
OR COUPLING UNIFIX 4 HOLE
REF. 1600363-001



TURBINES

TURBINES



BORA 2 LED
£ 485+VAT

2 YEARS WARRANTY

1 x BORA 2 LED
 REF. 1601152-001

- 20 Watts
- 59 dBa
- Head size - Ø 11,9 mm, h 13,2 mm
- 4 asymmetrical mixed air/water sprays
- LED Lighting



BORA 2 LK
£ 485+VAT

2 YEARS WARRANTY

1 x BORA 2 LK
 REF. 1601153-001

- 20 Watts
- 59 dBa
- Head size - Ø 11,9 mm, h 13,2 mm
- 4 asymmetrical mixed air/water sprays



BORALINA
£ 295+VAT

1 YEAR WARRANTY

1 x BORALINA
 REF. 1600373-001

- 12 Watts
- 62 dBa
- 1 spray (separate 2 x1)



BLACK PEARL ECO 4V
£ 170+VAT

1 YEAR WARRANTY

1 x BLACK PEARL ECO 4V
 REF. 1600040-001

- 13 Watts
- 62 dBa
- 1 spray (separate 2 x1)



DUO PACK BORA 2 LED
+ FREE COUPLING
 (CHOOSE BETWEEN UNIFIX OR MULTIFLEX)

£ 890+VAT

2 x BORA 2 LED
 REF. 1601152-001

+ COUPLING LK
 REF. 1600902-001
 OR COUPLING UNIFIX 4 HOLE
 REF. 1600363-001



DUO PACK BORA 2 LK
+ FREE COUPLING
 (CHOOSE BETWEEN UNIFIX OR MULTIFLEX)

£ 890+VAT

2 x BORA 2 LK
 REF. 1601153-001

+ COUPLING LK
 REF. 1600902-001
 OR COUPLING UNIFIX 4 HOLE
 REF. 1600363-001



COUPLINGS



COUPLING UNIFIX 4 HOLE (ELECTRIFIED)

£ 175+VAT

REF 1600363-001



COUPLING RA LK

£ 210+VAT

REF 1600902-001



COUPLING UNIFIX 4 HOLE (ELECTRIFIED)

£ 230+VAT

REF 1600086-001

- Spray Adjuster



COUPLING

£ 130+VAT

REF 1600082-001

- without light



COUPLING RA LK

£ 235+VAT

REF 1600866-001

- Spray Adjuster

NEW

NOVA

The Nova's small head and slim handle enhance your field of vision and guarantee easier access to the back of the mouth. The Nova has the **smallest head** and the **slimmest handle** in the entire Bien-Air handpiece range.

SEE BEYOND

MORE VISION, BETTER ACCESSIBILITY



CA NOVA 1:1 L REF 1601137-001
CA NOVA 1:1 L MICRO-SERIES REF 1601136-001

CA NOVA 1:5 L REF 1601139-001
CA NOVA 1:5 L MICRO-SERIES REF 1601138-001



Head 4x more resistant

Thanks to the new design, the Nova offers double the service life of the previous generation and makes it **4 times** more durable.



Extended service life

The **Cylcro gears** are designed to reduce friction and protect against wear. This new design makes the Nova an ultra-durable handpiece.



Improved reliability

The chuck system now features a new **PVD coating**. This makes the opening/closing system more reliable and extends the service life by 30 % compared to a conventional handpiece.




Enhanced safety

Integrated in the ceramic push-button thermally insulated from the Nova's head, the patented **Cool Touch +** heat-arresting technology prevents the build-up of heat to reduce the risk of burns to patients.



Say goodbye to obstacles

To guarantee optimal irrigation in the working area, the **Accu-Spray Quattro Mix** system focuses four asymmetrical air/water sprays* towards the end of the bur, with laser precision.



Anti-return valve

Risks of cross contamination are reduced thanks to the non-return valve and the **Sealed Head** protection. The non-return valve considerably reduces the flow of fluids from the oral cavity back into the dental unit's hoses



CA Classic

CA NOVA 1:5 L



AVAILABLE IN
micro
SERIES

NOVA 1:5L
£ 685^{+VAT}

1X NOVA 1:5L
REF 1601139-001
1X NOVA 1:5L MICRO-SERIES
REF 1601138-001



AVAILABLE IN
micro
SERIES

NOVA 1:1L
£ 545^{+VAT}

1X NOVA 1:1L
REF 1601137-001
1X NOVA 1:1L MICRO-SERIES
REF 1601136-001

CONTRA-ANGLES – HANDPIECES



AVAILABLE IN
micro
SERIES

CA 1:5L

£ 580+VAT
2 YEARS WARRANTY

1X CA 1:5L MICRO-SERIES
REF 1600690-001
OR 1X CA 1:5L STANDARD
REF 1600386-001



AVAILABLE IN
micro
SERIES

CA 1:1L

£ 465+VAT
2 YEARS WARRANTY

1X CA 1:1L MICRO-SERIES
REF 1600691-001
OR 1X CA 1:1L STANDARD
REF 1600384-001



AVAILABLE IN
micro
SERIES

PM 1:1

£ 265+VAT
2 YEARS WARRANTY

1X PM 1:1 MICRO-SERIES
REF 1600693-001
OR 1X PM 1:1 STANDARD
REF 1600383-001



CA ENDO 4:1

£ 630+VAT
1 YEAR WARRANTY

1X CA ENDO
REF 1600955-001



CAP PROPHY 15:1

£ 470+VAT
2 YEARS WARRANTY

1X CAP 15:1
REF 1600290-001



PM 10:1 PROPHY

£ 275+VAT
2 YEARS WARRANTY

1X PM 10:1 PROPHY
REF 1600289-001

* Disposable prophy heads
not sold with handpiece

MICROMOTORS



MICROMOTOR
MC3 LED / MC3 LK

£ 755+VAT
18 MONTHS WARRANTY

1X MC3 LED REF 1600680-001
OR 1X MC3 LK REF 1600077-001



MICROMOTOR
MC2 LED / MC2 ISOLITE 300

£ 730+VAT
18 MONTHS WARRANTY

1X MC2 LED REF 1600681-001
OR 1X MC2 ISOLITE 300 REF 1600078-001



MICROMOTOR
MX-i LED 3RD GEN

£ 805+VAT
3 YEAR WARRANTY

1X MX-i LED REF 1601008-001



MICROMOTOR
MX-i LED PLUS

£ 850+VAT
3 YEAR WARRANTY

1X MX-i LED REF 1600755-001

FOR IMPLANTOLOGY



CA 20:1 L

£ 985+VAT

2 YEARS WARRANTY

1 x CA 20:1L MICRO-SERIES
REF. 1600692-001
OR 1 x CA 20:1L
REF. 1600598-001

AVAILABLE IN
micro
SERIES



PM 1:1

£ 295+VAT

2 YEARS WARRANTY

1 x PM 1:1
REF. 1600052-001

- External spray
- Compatible with all the range of micromotors



SET MX-i LED PLUS + CABLE
(FOR ICHIROPRO AND CHIROPRO PLUS)

£ 1 182+VAT ~~£ 1 295+vat~~

REF 1700758-001

-50%
on the cable



SET MX-i LED + CABLE
(FOR CHIROPRO 3RD GEN)

£ 1 120+VAT ~~£ 1 250+vat~~

REF 1700744-001

-50%
on the cable

OFFER

MAINTENANCE PACK

FAST, EFFICIENT AND ECONOMIC MAINTENANCE

Lubricare 2 performs handpiece maintenance (lubrication and cleaning), which is essential for the operation of instruments. Up to 4 handpieces can be serviced automatically in a very short time. The amount of lubricant and the maintenance time for each step are pre-programmed in the Lubricare 2

>>> **Up to 4 instruments at once in no time at all**


 **ONLY 20 SECS.**

Thanks to the different pre-programmed operating modes, adapted specifically to each handpiece, Lubricare 2 can automatically perform maintenance of four instruments in a short time. Lubricare 2 effectively removes excess oil, which allows for quick maintenance. It cleans a turbine in about 20 seconds and a contra-angle in 40 seconds

>>> **Easy cleaning**

The control panel is smooth to facilitate cleaning and the door can be removed for a better access. Lubricare 2 incorporates a mechanism to extract the oil vapors and prevent the diffusion of aerosols during its use

>>> **Spraynet® your best friend**

 **99.99% reduction of the microbial load**

To obtain better results, Bien-Air strongly recommends the use of Spraynet®, proven to reduce 99.99% of the microbial load.

>>> **Economic maintenance**

4X less oil consumption

When using Lubricare 2, oil consumption is reduced by 4 times compared to manual procedure, thus allowing for a faster maintenance with excellent cost savings.

LUBRICARE 2

£ 1 800 +VAT

1 X LUBRICARE 2 REF. 1601156-001
1 X DUO PACK MAINTENANCE
LUBRIFLUID REF 1600064-001
SPRAYNET REF 1600036-001



SWISS  MADE

FOR AN ADDITIONAL **£ 180 +VAT.**

FOR **£ 1 980 +VAT**

RECEIVE

2 COUPLINGS OF YOUR CHOICE

TURBINE COUPLINGS FOR LUBRICARE 2 :

- W&H
- Sirona
- KaVo
- Bien Air



CONSUMABLES



LUBRIFLUID
BOX OF 6 BOTTLES
£ 138 +VAT
REF 1600064-006



SPRAYNET
BOX OF 6 BOTTLES
£ 138 +VAT
REF 1600036-006



AQUACARE
BOX OF 6 BOTTLES
£ 95 +VAT
REF 1600617-006



LUBRIMED CATRIDGES
£ 19 +VAT
REF 1600037-006



BOX OF 10 IRRIGATION LINES
£ 66 +VAT
REF 1500984-010

CONTACT YOUR DEALER OR BIEN-AIR UK FOR CONSUMABLES



iCHIROPRO

ULTIMATE IMPLANTOLOGY

IMPLANTOLOGY

ORAL SURGERY

The unique and revolutionary iChiropro is the perfect partner for your implantology and dental surgery operations.

Controlled by an iPad application, this intuitive system supports you in every step of your procedures. Simply select the pre-programmed operation of your choice and the iChiropro will provide instructions and guidance. For optimal patient file analysis and management, the application can be used to import data from the coDiagnostiXTM software and you can also save and export your operation reports with a single click.

Simplicity, comfort and modernity: discover a new work philosophy!

S W I S S  M A D E

SET iCHIROPRO

FROM

£ 3 275_{+VAT}

1 x SET iCHIROPRO REF. 1700439-001

SET iCHIROPRO with 1 contra-angle CA 20:1 L

£ 3 960_{+VAT}

1 x SET iCHIROPRO REF. 1700440-001

SET iCHIROPRO with 2 contra-angles CA 20:1 L

£ 4 860_{+VAT}

1 x SET iCHIROPRO REF. 1700440-001
+ 1 X CA 20:1 L MICRO-SERIES REF. 1600692-001

The offer includes :

- 1 x console iChiropro
- 1 x micromotor MX-i LED
- 1 x cable MX-i LED
- 1 x pedal
- 1 x handpiece stand
- 10 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

The offer includes :

- 1 x console iChiropro
- 1 x contra-angle CA 20:1L
- 1 x micromotor MX-i LED
- 1 x cable MX-i LED
- 1 x pedal
- 1 x handpiece stand
- 10 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

The offer includes :

- 1 x console iChiropro
- 2 x contra-angles CA 20:1L
- 1 x micromotor MX-i LED
- 1 x cable MX-i LED
- 1 x pedal
- 1 x handpiece stand
- 10 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY



iPad not included

- ✓ Multi-implant mode
- ✓ Pre-programmed operation sequences
- ✓ Operation report and patient file
- ✓ Linked to CoDiagnostix (implantology planning software)
- ✓ Multiple Users
- ✓ Implant barcode reader
- ✓ Torque > 70 Ncm

CHIROPRO PLUS

IMPLANTOLOGY

ORAL SURGERY



The Chiropro Plus has been especially developed for oral surgeons and maxillofacial surgeons. Thanks to the many possible motor-instrument combinations, the Chiropro PLUS is the perfect solution for implant placement, periodontology and oral surgery, including wisdom tooth extraction.



SWISS  MADE

CHIROPRO PLUS

FROM

£ 2 575_{+VAT}

1 x SET CHIROPRO PLUS 3rd GEN REF. 1700710-001

The offer includes :

- 1 x console Chiropro PLUS 3rd Gen
- 1 x micromotor MX-i LED PLUS 3rd Gen with 3 year warranty
- 1 x cable MX-i LED PLUS 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

CHIROPRO PLUS with 1 contra-angle CA 20:1 L

£ 3 190_{+VAT}

1 x SET CHIROPRO PLUS 3rd GEN REF. 1700709-001

The offer includes :

- 1 x console Chiropro PLUS 3rd Gen
- 1 x contra-angle CA 20:1 L Micro-Series
- 1 x micromotor MX-i LED PLUS 3rd Gen with 3 year warranty
- 1 x cable MX-i LED PLUS 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

CHIROPRO PLUS with 2 contra-angles CA 20:1 L

£ 4 090_{+VAT}

1 x SET CHIROPRO PLUS 3rd GEN REF. 1700709-001
+ 1 X CA 20:1 L MICRO-SERIES REF. 1600692-001

The offer includes :

- 1 x console Chiropro PLUS 3rd Gen
- 2 x contra-angles CA 20:1 L Micro-Series
- 1 x micromotor MX-i LED PLUS 3rd Gen with 3 year warranty
- 1 x cable MX-i LED PLUS 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

✓ Clear and concise interface

Sober and simple, at a glance the screen gives you access, to all the information essential to the smooth running of your procedure: type of instruments, speed, torque, irrigation flow and direction of rotation

✓ Easy to maintain

The control button can be easily removed from the console and sterilized. Moreover, the screen is protected by glass, highly resistant and easy to clean.

✓ Intuitive Navigation

Thanks to its unique control button, navigation becomes child's play. Turn the button to scroll through the different menus and then validate your choice with a simple press. You said simplicity?

✓ Absolute precision

Thanks to an accurate high-tech electronic system, the torque selected on the console corresponds exactly to the torque obtained at the output of your instrument. No additional calibration is therefore required. It's as simple as it appears.

✓ Preprogrammed operating sequences.

Thanks to the preprogrammed operating sequences, you can carry out your procedures in a minimum of time. Depending on your patient's clinical data, determine the number of steps necessary for the realization of your procedure.

CHIROPRO

General dental practitioner or experienced surgeon. Chiropro is the simple and fast solution for all your implant procedures.

Thanks to its unique control button, navigation is easier than ever before. Turn the button to scroll through the different menus and then validate your choice with a simple press.

You said simplicity?

IMPLANTOLOGY

 PURE SIMPLICITY



SWISS  MADE

CHIROPRO

FROM

£ 2 165 +VAT

1 x SET CHIROPRO 3RD GEN REF. 1700708-001

The offer includes :

- 1 x console Chiropro
- 1 x micromotor MX-i LED 3rd Gen
- 1 x cable MX-i LED 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

CHIROPRO

with 1 contra-angle
CA 20:1 L

£ 2 750 +VAT

1 x SET CHIROPRO 3RD GEN REF. 1700707-001

The offer includes :

- 1 x console Chiropro
- 1 x contra-angle CA 20:1 L Micro-Series
- 1 x micromotor MX-i LED 3rd Gen
- 1 x cable MX-i LED 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

CHIROPRO

with 2 contra-angles
CA 20:1 L

£ 3 650 +VAT

1 x SET CHIROPRO 3RD GEN REF. 1700707-001
+ 1 X CA 20:1 L MICRO-SERIES REF. 1600692-001

The offer includes :

- 1 x console Chiropro
- 2 x contra-angles CA 20:1 L Micro-Series
- 1 x micromotor MX-i LED 3rd Gen
- 1 x cable MX-i LED 3rd Gen
- 1 x pedal
- 1 x handpiece stand
- 5 x irrigation lines
- 1 x trio maintenance pack

2 YEARS WARRANTY

✓ Driven by **Chiropro**, the combination of the new MX-i micromotor and CA 20:1 contra-angle handpiece offer you the best of rotary technology for all your implantology procedures.

✓ Most popular and durable implantology handpiece on the market*, the CA 20:1 contra-angle has an internal irrigation system. So you are not bothered by the irrigation line when holding the instrument.

✓ Coupled with the MX-i micromotor, the CA 20:1 contra-angle handpiece (also available in Kirschner/Meyer and Micro-Series versions) ensures unparalleled speed stability, for smooth and precise interventions without vibration.

✓ Made easier, your procedures can be performed with confidence.

GO ELECTRIC WITH BIEN-AIR INTEGRATED MICROMOTOR KITS



SET iOPTIMA^{INT}

£ 2 850^{+VAT}

1 x SET iOPTIMA INT REF 1700704-001
+ REF 6200217-001 (INCLUDES IPAD MINI)

- All-in-One solution Offering full range of treatment. Restorative, endodontics (Rotary Endo), piezo scaling and implantology*
- WaveOne Gold and Reciproc (Optional Purchase of additional H/P Required)
- Numerous pre-set and customizable protocols including Non-AGP
- Quickly integrated into almost any existing dental unit.

* Controls piezo scalers, scaler not included. Implant motor and pump optional extra, not included.

MX2 OPTIMA INT PRO

£ 2 150^{+VAT}

1 x MX2 OPTIMA INT PRO REF 6200186-001

£ 1 955^{+VAT}

1 x MX2 OPTIMA INT PRO W/O TRANSFORMER REF 6200164-001



- Programmable console
- Non-AGP Mode
- Full endo motor (rotary)
- Led Light
- 3 Years Warranty on MX2 Micromotor
- Semi-Swivel Hose
- Reverse
- Power/Torque - 3.5Ncm

MCX OPTIMA INT

£ 1 675^{+VAT}

1 x MCX OPTIMA INT REF 6200222-001

£ 1 425^{+VAT}

1 x MCX OPTIMA INT W/O TRANSFORMER REF 6200183-001



- Programmable console
- Non-AGP Mode
- Led Light
- 3 Years Warranty on MCX Motor
- Non-Swivel Hose
- Reverse
- Power/Torque - 2.5Ncm

MCX LED

£ 1 170^{+VAT}

1 x MX2 LED REF 6200188-001

£ 995^{+VAT}

1 X MCX LED W/O TRANSFORMER REF. 6200189-001



- Control thru foot pedal
- Speed limited by Potentiometer
- LED Light
- 3 Years Warranty on MCX Micromotor
- Non-Swivel Hose
- No Reverse
- Power/Torque - 2.5Ncm

OPTIMA MCX

(SEMI-INTEGRATED MICROMOTOR SOLUTION)

£ 1 215^{+VAT}

1 x OPTIMA MCX REF. 1700595-001
1 x OPTIMA MCX BLUE REF. 1700596-001
1 x OPTIMA MCX ORANGE REF. 1700597-001
1 x OPTIMA MCX GREEN REF. 1700598-001
1 x OPTIMA MCX PINK REF. 1700599-001



- Easily connects to most dental chairs
- Control thru foot pedal/Set dial for speed
- LED Light
- 3 Years Warranty on MCX Micromotor
- Non-Swivel Hose
- Power/Torque - 2.5Ncm