

# Summer Offers



VALID UNTIL  
31/08/2023



Endo



Surgical



Decon



Restorative



Oral Hygiene

**Support & Service**

**and  
much  
more!**

## NLX nano RETRO FITTED MICROMOTOR

Electric micromotors are the driving force behind speed increasing handpieces and can be fully integrated into the delivery system of most new and existing dental chairs. They can be retro fitted by an engineer or come as an option when buying new chairs. Controlling the specific speed and torque of micromotors within existing platforms is made easy thanks to built-in control pads.

Speed increasing contra-angles not only do everything a turbine does using the same FG burs, they provide a constant torque from the micromotor that does not decrease with resistance and higher load with no loss of power as they come in contact with the tooth.


Scan/Click

For more info




To discuss the offers & options contact NSK:

**OFFER 1 BUY 3 Ti-Max Z95L AND GET OPTION 1 FREE**



**OFFER 2 BUY 4 Ti-Max Z95L AND GET OPTION 3 FREE**



**OFFER 3 BUY ANY OF THE BUILT-IN OPTIONS BELOW AND GET FREE Ti-Max Z95L SPEED INCREASING CONTRA-ANGLE**



- option 1 NLX nano micromotor & BUILT-IN KIT \*** (including MultiPad) RRP £2,229
- option 2 NLX nano micromotor & BUILT-IN KIT \* + Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling perio, endo & restorative) RRP £3,531
- option 3 NLX plus micromotor** (includes Rotary Endo functions) & BUILT-IN KIT \* (including MultiPad) RRP £2,681
- option 4 NLX plus micromotor** (includes Rotary Endo functions) & BUILT-IN KIT \* (including MultiPad) + **Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling, perio, endo & restorative) RRP £3,957
- option 5 Varios170 Lux** ULTRASONIC HANDPIECE (tip options for scaling perio, endo & restorative) + BUILT-IN KIT \* (including MultiPad) RRP £1,495

\*Excludes install. Through participating dealers only.

## NLX nano PORTABLE MICROMOTOR

- Portable electric micromotor for connecting slow, straight and speed increasing handpieces
- Puts you in control - set your own speeds
- Helps to manage aerosol generation
- Connects to your existing cart

MODEL NLX nano Portable  
CODE Y1001362  
RRP £2,629



**OFFER NLX nano Portable ONLY £1,656 + FREE Ti-Max Z95L, WORTH £1,556**



# Contra-angles

## Ti-Max

2 + 1 FREE\*

- Lightweight external titanium body • **Optic with internal water** • Small head/narrow neck
- **Unique spray switch on Z95L, Z45L, Z25L** • 3 year warranty



- **Optic**
- **1:5 Increasing+Spray switch**
- Head size  $\varnothing 8.9 \times 13.4\text{mm}$

MODEL Z95L  
CODE C1034001-SW  
RRP £1,556



- **Optic**
- **1:5 Increasing**
- 10% shorter & lighter than Z95L
- Head size  $\varnothing 8.9 \times 13.4\text{mm}$

MODEL nano95LS  
CODE C1099001  
RRP £1,482



- **Optic**
- **1:4.2 Increasing+Spray switch**
- 45° head
- Head size  $\varnothing 9.6 \times 14.5\text{mm}$
- Ideal for minor oral surgery too
- For regular or extra long (25mm) FG burs

MODEL Z45L  
CODE C1064001  
RRP £1,556



- **Optic**
- **1:1 Direct Drive +Spray switch**

MODEL Z25L  
CODE C1038001-SW  
RRP £1,233

Scan/Click  
For a full range of  
NSK Contra-angles

## S-Max

3 + 1 FREE\*

- Stainless steel • **All with internal water** • Choice of optic & non-optic • 2 year warranty



- **Optic**
- **1:5 Increasing**
- Head size  $10 \times \text{H.}14.6\text{mm}$

MODEL M95L  
CODE C1023001  
RRP £1,005

- **Non-Optic**
- **1:5 Increasing**
- Head size  $10 \times \text{H.}14.6\text{mm}$

MODEL M95  
CODE C1026001  
RRP £834



- **Optic**
- **1:1 Direct Drive**

MODEL M25L  
CODE C1024001  
RRP £814

- **Non-Optic**
- **1:1 Direct Drive**

MODEL M25  
CODE C1027001  
RRP £535



- **Optic**
- **4:1 Reducing**

MODEL M15L  
CODE C1025001  
RRP £952

- **Non-Optic**
- **4:1 Reducing**

MODEL M15  
CODE C1028001  
RRP £636



- **Non-Optic**
- **1:1 Straight**

MODEL M65  
CODE H1008001  
RRP £523

Scan/Click  
For a full range of  
NSK Contra-angles

## FX

1 + 1 FREE\*

- 12 month warranty



- **Non-Optic Contra-angle**
- **1:1 Direct Drive**
- **No internal water**
- Push button

MODEL FX25m  
CODE C1052001  
RRP £264



- **Non-Optic Straight**
- **1:1 Direct Drive**
- **No internal water**

MODEL FX65m  
CODE H1014001  
RRP £264

\*Terms and conditions apply. Offers only available through participating dealers.

## NAC-EC

3 FOR £195\*

- 12 month warranty



- **Non-Optic Contra-angle**
- **1:1 Direct Drive**
- **No internal water**

MODEL NAC-EC  
CODE Y110148  
RRP £155

Scan/Click  
For a full range of  
NSK Contra-angles

## Ti-Max

2 + 1 FREE\*

• Lightweight external titanium body • Clean head system • 4 hole spray • **3 year warranty**

NEW



- **Optic**
- **Power 26w**
- **Std head:**  
Ø 12.5 x 13.1 mm

MODEL Z900L  
CODE P1111001  
RRP £1,134

Also available in fittings for  
KaVo®, Sirona®, W&H®, Bien-Air®  
RRP £1,191



- **Optic**
- **Power 23w**
- **Mini head:**  
Ø 10.8 x 12.1 mm

MODEL Z800L  
CODE P1110001  
RRP £1,134

Also available in fittings for  
KaVo®, Sirona®, W&H®, Bien-Air®  
RRP £1,191



- **Optic**
- **Power 20w**
- **Micro head:**  
Ø 9 x 10.8 mm

MODEL Z Micro  
CODE PA2357001  
RRP £1,134

Also available in fittings for  
KaVo®, Sirona®  
RRP £1,191  
for regular or short shank burs

Scan/Click  
For full info on  
models/order codes

## S-Max

OPTIC 3 + 1 FREE\*

NON-OPTIC 3 FOR £1,122\*

• Stainless steel • Enhanced power • 4 hole spray • DIY cartridge change • Clean head system • **2 year warranty**



- **Optic**
- **Power 26w**
- **Std head:**  
Ø 12.1 x 13.3 mm

MODEL M900L  
CODE P1254001  
RRP £747

Also available in fittings for  
KaVo®, Sirona®, W&H®, Bien-Air®



- **Optic**
- **Power 23w**
- **Mini head:**  
Ø 10.6x12.4mm

MODEL M800L  
CODE P1253001  
RRP £747

Also available in fittings  
for KaVo®



- **Non-Optic**
- **Power 26w**
- **Std head:**  
Ø 12.1 x 13.3 mm

MODEL M900  
CODE P1256001  
RRP £696

Also available in fittings  
for KaVo®, W&H®



- **Non-Optic**
- **Power 23w**
- **Mini head:**  
Ø 10.6x12.4mm

MODEL M800  
CODE P1255001  
RRP £696

Scan/Click  
For full info on  
models/order codes

## Pana-Max2

1 + 1 FREE\*

• Non-Optic turbine with single hole spray • Save time/money with DIY cartridge change • Clean head system • **12 month warranty**



- **Non-Optic**
- **Head:**  
Ø 11.2 x 13.4 mm
- **Coupling connection**

MODEL Pana-Max2  
CODE P1211001  
RRP £453

for NSK PTL coupling

- **Non-Optic**
- **Head:**  
Ø 11.2 x 13.4 mm
- **Coupling connection**

MODEL Pana-Max2 KV  
CODE P1212001  
RRP £487

for KaVo® MULTIflex® coupling



- **Non-Optic**
- **Head:**  
Ø 11.2 x 13.4 mm
- **Direct connection**

MODEL Pana-Max2 M4  
CODE P1154001  
RRP £411

For Midwest M4 hose

- **Non-Optic**
- **Head:**  
Ø 11.2 x 13.4 mm
- **Direct connection**

MODEL Pana-Max2 B2  
CODE P1155001  
RRP £411

for Borden hose

Scan/Click  
For a full list of  
fittings/models/codes

\*Terms and conditions apply. Offers only available through participating dealers.

## VarioSurg3

ULTRASONIC BONE SURGERY SYSTEM

**50% OFF** FOR EXISTING  
**Surgic Pro / Pro2 / XT USERS\*†**

OR

**FREE 2<sup>nd</sup> HANDPIECE/CORD**  
**OR 12 EXTRA TIPS**  
**OR AN OSSEO 100+**

- Piezo technology bone surgery system with LED.
- Lightweight slimline handpiece.
- Contains a tip stand with 6 tips – SG11, SG6D & SG7D for sinus lifts.  
H-SG1 for cutting bone with SG3 & SG5 for scraping.
- Wide range of tips available for fast and more precise treatments.

MODEL VarioSurg3  
CODE Y1002726  
RRP £8,086

Scan/Click  
For a full range of tips

Scan/Click  
For the core range of tips



**VarioSurg tips 1 + 1 FREE\***



## LINK OPTIONS

### Link VarioSurg3 with Surgic Pro or Surgic Pro2

Control two units with one foot pedal. Choose a Link Stand or store everything on an iCart where you can adjust the height and move it around on its wheels.



#### Linking Kit for Surgic Pro

MODEL Link Stand/Cable  
CODE Y1002729  
RRP £512



#### Linking Kit for Surgic Pro2

MODEL Link Stand2  
CODE Z1452001  
RRP £301

MODEL Wireless Link Module  
CODE Z1402001  
RRP £237



#### iCart

MODEL iCart Duo  
CODE S9090  
with Surgic Pro link cable (CODE Z1189001)

MODEL iCart Duo  
CODE S9090  
with Surgic Pro2 link module (CODE Z1402001)

RRP £723

## Surgic Pro2

**NEW**

### SURGICAL MICROMOTOR

- Compact new motor: 25% shorter, 9% lighter than current motor – Greatly improved balance
- Link by Bluetooth
- Significantly improved optics (LED) • **2 year warranty**

Complete Kit contains:

- Control unit
- Motor/Cord
- 20:1 contra angle
- Wireless foot control

- **Optic**
- **With Ti-Max X-SG20L**

MODEL Surgic Pro2  
CODE Y1004195  
RRP £5,856

- **Non-Optic**
- **With S-Max SG20**

MODEL Surgic Pro2  
CODE Y1004196  
RRP £4,164



**BUY ANY Surgic Pro2 AND GET Osseo 100+ FOR FREE\*\***



MODEL Osseo 100+  
CODE Y1004176

## LINK BY BLUETOOTH

- Wireless Foot Control
- iPad
- Variosurg3 (Surgical Piezo) operate with 1 wireless Foot control
- Osseo100+ confirms ISQ value over a shared interface



## ADDITIONAL MOTOR/CORD



- **Optic**

MODEL SGL80M  
CODE Y1004211  
RRP £2,476



- **Non-Optic**

MODEL SG80M  
CODE Y1004212  
RRP £1,772

**FREE 20:1 IMPLANT HANDPIECE OR SURGICAL STRAIGHT**  
WHEN YOU BUY AN ADDITIONAL MOTOR/CORD



## Osseo 100+

### OSSEOINTEGRATION MONITORING DEVICE

Osseo 100+ measures implant stability and osseointegration to enhance decisions about when to load the implant. This is invaluable when using protocols with shorter treatment times and treating higher risk patients. The simple steps to measure the ISQ value allow the implant loading period to be planned in advance.

MODEL Osseo 100+  
CODE Y1004176  
RRP £2,570

**£1,799**



Scan/Click  
More info &  
Multipeg List

\*Terms and conditions apply. Offers only available through participating dealers.

## Surgic Pro

### SURGICAL MOTOR UNIT

Complete Kit contains:

- Control unit
- Motor/Cord
- 20:1 contra angle
- **2 year warranty**
- Corded foot control

- **Optic**
- **With Ti-Max X-SG20L**

MODEL Surgic Pro  
CODE Y1003586  
RRP £6,002

- **Non-Optic**
- **With S-Max SG20**

MODEL Surgic Pro  
CODE Y1003587  
RRP £4,481



**FREE 2 ADDITIONAL 20:1 CONTRA-ANGLES\***  
WITH Surgic Pro COMPLETE KITS

### ADDITIONAL MOTOR/CORD



- **Optic**

MODEL SGL70M  
CODE E1023001  
RRP £2,665



- **Non-Optic**

MODEL SG70M  
CODE E1025001  
RRP £1,853

### IRRIGATION TUBING

- Used with Surgic Pro / Pro2 / XT and VarioSurg3



Irrigation tubing x5 (disposable)  
CODE Y900079  
RRP £62.20



Irrigation Tubing x1 (Autoclavable)  
CODE S900055  
RRP £52

**Eco Friendly  
autoclavable  
option**

**FREE 20:1 IMPLANT HANDPIECE  
OR SURGICAL STRAIGHT**  
WHEN YOU BUY A SPARE MOTOR/CORD




### ZYGOMATIC HANDPIECE

**NEW**

- **Non-Optic** • **20:1 Reduction** • **Twist chuck** • **12 Month warranty**

NSK's new handpiece for Zygomatic & Pterygoid implants. Its unique angulated shape offers optimal access and visibility.

MODEL SGX-E20R  
CODE HA1200001  
RRP £1,529



**1 + 1 FREE\***

### IMPLANT CONTRA-ANGLES 20:1

- Max speed: 2,000min<sup>-1</sup> • **2 year warranty**



- **Optic & Dismantleable**

MODEL Ti-Max X-DSG20L CODE C1068001 RRP £1,527

- **Optic**

MODEL Ti-Max X-SG20L CODE C1003001 RRP £1,243

- **Non-Optic**

MODEL S-Max SG20 CODE C1010001 RRP £1,049

### SURGICAL STRAIGHTS 1:1

- Used with standard HP burs • **2 year warranty**



- **Optic**

MODEL Ti-Max X-SG65L CODE H1009001 RRP £933

**TRIPLE PACK OFFER £1,855\***

- **Non-Optic**

MODEL Ti-Max X-SG65 CODE H1038001 RRP £636

**TRIPLE PACK OFFER £1,308\***

\*Terms and conditions apply. Offers only available through participating dealers.

# Mobile dentistry

## VIVA ace

### MOBILE DENTISTRY SYSTEM

- Ideal as a surgery back up or to visit shielded patients
- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic handpiece
- Quick setup and storage
- 4-level adjustable height and angle



MODEL VIVA ace  
CODE Y1003773  
RRP £12,517

Weight 8.6kg

Width 462 mm, Height 284 mm, Depth 292 mm

## iProphy

NEW

### CORDLESS PROPHY HANDPIECE

The complete kit comes with a 4:1 reduction head for screw-in prophy cups/brushes. Alternatively a PR-F head (optional extra) can be used that connects to disposable Prophy heads.



MODEL iProphy Complete Kit  
CODE Y1500015  
RRP £1,085

**FREE S-Max M95L, M25L & M65 Handpieces\***

**OFFER 2 ×**  **FREE**

Monthly **RENTAL** £433\* Fully maintained for 3 years

# Endodontics

## Endo-Mate TC2

### CORDLESS ENDO MOTOR/HANDPIECE

- Large LCD display
- 5 programmes for different file systems
- Supports most major brands of Ni-Ti files
- Torque control facility ranging from 0.1 to 4.0 Ncm depending on head:gear ratio
- Auto-reverse and alarm function sounds an audible alert when preset torque level has been reached

MODEL Endo-Mate TC2  
CODE Y1001027  
RRP £1,710



## Endo-Mate DT2

### CORDED ENDO MOTOR/HANDPIECE

- Large LCD display
- Store up to 9 speed and torque settings
- Use with Ni-Ti files from most major suppliers
- Torque control from 0.1 to 6.5 Ncm depending on head:gear ratio
- Auto-reverse on, auto-stop, autoreverse off

MODEL Endo-Mate DT2  
CODE Y1004230  
RRP £1,616



## iPex II

### APEX LOCATOR

- Accurate measurement of any root canal using SmartLogic technology
- Operates in all conditions, including dry, wet and bloody canals

MODEL iPex II  
CODE Y1002208  
RRP £1,015



**BUY Endo-Mate DT2 AND GET iPex II FREE WORTH £1,015\***

\*Terms and conditions apply. Offers only available through participating dealers.



## Varios Combi Pro iPiezo engine®

TWO ORAL CARE SYSTEMS IN ONE!



Combines ultrasonic treatments and powder therapy in one unit. Comes with a detachable Prophy chamber. Detachable Perio chamber available as an optional extra.

MODEL Varios Combi Pro  
CODE Y1002843  
RRP £4,977



Portable cart for Varios Combi Pro

MODEL iCart Prophy  
CODE S9095  
RRP £719

Scan/Click

**FREE Perio Set WORTH £1,190\***

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

For more info

## Varios 970 iPiezo engine®

ULTRASONIC SCALER

Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/ perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Units come with:

- 1 handpiece
- G4, G8 & G16 tips
- 3 tip wrenches
- sterilisation case

• **Optic**

MODEL Varios 970 LUX  
CODE Y1004161  
RRP £2,898

• **Non-Optic**

MODEL Varios 970  
CODE Y1004160  
RRP £2,494



**OFFER + 2 FREE**



**HANDPIECES & STERILISATION CASES\*\***

\*\*EXCLUDING EXTRA TIPS/TIP WRENCHES

Varios2 handpiece

• **Optic**

MODEL VA2-LUX-HP  
CODE E351050  
RRP £531



Varios2 handpiece

• **Non-Optic**

MODEL VA2-HP  
CODE E350050  
RRP £279



**TRIPLE PACK OFFER £1,077\***

**TRIPLE PACK OFFER £630\***

**SCALER TIPS 1 + 1 FREE\***



Scan/Click  
Full range of tips

Scan/Click  
To order

\*Terms and conditions apply. Offers only available through participating dealers.

## Scaler Tips


1 + 1 FREE\*

TIPS FOR VARIOS  
ULTRASONIC SCALER

Scan/Click  
Full range  
of tips

TIPS FOR AIR SCALERS

Scan/Click  
Tips for  
Air Scaler



## AIR SCALERS



Direct connection  
air scaler

- Oscillation frequency 5,600-5,850Hz
- Kit includes 3 scaler tips

MODEL AS2000  
CODE T698001 (Midwest M4 hose)  
CODE T697001 (Borden B2/B3 hose)  
RRP £477 each



Optic air scaler with  
coupling connection

- Oscillation frequency 6,200-5,850Hz
- Scratch resistant with titanium body
- Kit includes 3 scaler tips

MODEL Ti-Max S970L  
CODE T1005001\*\*  
RRP £1,339

\*\*also available in KaVo® fitting T1015001  
and Sirona® fitting T1013001, RRP £1,372

FREE NAC-EC WITH AS2000

FREE 2 x NAC-EC WITH S970

## PROPHY HANDPIECES

FX 1 + 1 FREE\*




- Max speed 1,000 rpm
- Brushes/Cups: **Screw in**

MODEL FX57m  
CODE C1054001  
RRP £297

- Max speed 40,000 rpm
- Brushes/Cups: **Disposable Prophy heads**

MODEL FX65m  
CODE H1014001  
RRP £264



Paste defence system

NAC-EC 3 FOR £195\*



AR-EC(k) AR-EC(s) NAC-EC

- Max speed 5,000 rpm
- Brushes/Cups: **Snap on cups**

MODEL AR-EC(k)  
CODE Y110165  
RRP £153

- Max speed 5,000 rpm
- Brushes/Cups: **Screw in**

MODEL AR-EC(s)  
CODE Y110164  
RRP £153

- Max speed 30,000 rpm
- Brushes/Cups: **Latch grip**

MODEL NAC-EC  
CODE Y110148  
RRP £155

\*Terms and conditions apply. Offers only available through participating dealers.

## Prophy-Mate neo SUPRAGINGIVAL AIR POLISHER



- Compact and lightweight
- Includes 60° and 80° nozzles
- 360° rotation
- Packed in a handy case
- Easily dismantled for autoclaving
- Includes 6 sachets of FLASH pearl® that quickly removes stains and plaque

CODE Y135029 for NSK PTL coupling\*  
 CODE Y135030 for KaVo® MULTIflex® LUX coupling  
 CODE Y135031 for Sirona® quick coupling  
 CODE Y135032 for W&H® Roto Quick® coupling  
 CODE Y135033 for Bien-Air® Unifix® coupling

RRP £1,124



+  **FREE SECOND HANDPIECE WORTH £427\***



### Prophy-Mate neo nozzles

CODE T797012 60° nozzle  
RRP £120

CODE T797013 80° nozzle  
RRP £120

### FLASHpearl

- Calcium carbonate (54µm) based powder for stain removal



4 x 300 g bottles  
CODE Y900698  
RRP £138

100 x 15 g sachets  
CODE Y900693  
RRP £189

### SOFTpearl **NEW**

- Glycine based powder (65µm) for moderate stain removal. Pleasant taste with no flavourings, ideal for use on children, and patients with allergies.



4 x 250 g bottles  
CODE Y1500721  
RRP £133

**1 + 1 FREE\***

**BUY ANY 2 PACKS GET A Prophy-Mate neo NOZZLE FREE\***

## Perio-Mate BIOFILM ERASER



- Avoids injury to delicate periodontal pockets
- Smooth 360° rotation
- Flexible plastic nozzle tip for better access to difficult to reach areas
- Ultra-fine water spray and Perio-Mate powder eradicate biofilm in just 20 seconds
- Packed in a handy case with 160g powder

CODE Y1002654 for NSK coupling  
 CODE Y1002655 for KaVo® MULTIflex® LUX coupling  
 CODE Y1002656 for Sirona® quick coupling  
 CODE Y1002657 for W&H® Roto Quick® coupling  
 CODE Y1002658 for Bien-Air® Unifix® coupling

RRP £1,281



+  **FREE SECOND HANDPIECE WORTH £427\***



- Glycine based powder (25µm)
  - 2x 160 g bottles
- Perio-Mate Powder  
CODE Y900938  
RRP £96.10



- Pack of 40 Perio-Mate Tip Nozzles  
CODE Y1002741  
RRP £96

**Prophy-Mate neo + Perio-Mate**

**£1,854\***

# SMART<sup>3</sup> Decontamination Solutions


AFFORDABLE,  
RELIABLE  
& RAPID VACUUM  
STERILISATION



## COPPER CHAMBER

Save on running costs

- Copper chamber more efficient to run than traditional stainless steel
- Quicker instrument turnaround
- Effective drying – protects from corrosion and recontamination
- Improved reliability – reduced maintenance costs
- Extends life of equipment and handpieces



Scan/Click

For more info

**EXTENDED WARRANTY  
ON NSK HANDPIECES  
WITH EVERY NEW NSK  
AUTOCLAVE**

## iCare



Automated handpiece lubrication unit.

Compressed air connection required 5-6 bar.  
6mm air line required, ideally with a shut off valve.

MODEL C2  
CODE Y1002787  
RRP £2,570

MODEL C3  
CODE Y1002788  
RRP £2,701

**HALF PRICE £1,285**  
or **£1,065** when buying a **NSK autoclave**

**RENT** for just £26.50pm<sup>†</sup>

## iMax S



Rapid, reliable S-type autoclave.

- 20 min fast cycles for unwrapped loads
- 25 min handpiece cycle
- Affordable, effective upgrade from non-vacuum autoclaves
- Full vacuum drying phase (Dry Handpieces!)

CODE DX91803S0200  
RRP £3,070

**ONLY £3,070**  
**£350\*\* install & delivery**

**RENT** for just £115pm<sup>†</sup>

## iClave plus



High performance B-type vacuum autoclave.

- Easy to use
- Fast 25 min cycle for unwrapped loads
- 35 min for pouched loads
- Ideal for implant kits
- Includes Helix Test + Strips

CODE DXY1003077  
RRP £6,197

**FROM ONLY £3,995**  
**£350\*\* install & delivery**

**RENT** for just £155pm<sup>†</sup>

<sup>†</sup>Visit [mynsk.co.uk/rental](http://mynsk.co.uk/rental) for more information. <sup>\*\*</sup>For Scotland, call for details.

## PROTECT IT

All-Inclusive Maintenance and Breakdown Contracts for NSK autoclaves from just **£64pm.**

One annual maintenance visit or every 1500 cycles	Additional maintenance visit for units exceeding 1500 cycles	All service parts and service sundries included	All breakdown calls*	All breakdown parts including labour*	All engineer travel in the UK	Free on-site operator training	Free technical software upgrades
---	--	---	----------------------	---------------------------------------	-------------------------------	--------------------------------	----------------------------------

## CONTACT NSK

**Scotland - Sales Manager**  
Angela Glasgow  
[glasgow@nsk-uk.com](mailto:glasgow@nsk-uk.com) 07525 911006

**Ireland - Sales & Services Manager**  
Jonathan Singh  
[jsingh@nsk-uk.com](mailto:jsingh@nsk-uk.com) 07464 675158

**Public Sector & Decon Manager**  
Mark Beckwith  
[mbeckwith@nsk-uk.com](mailto:mbeckwith@nsk-uk.com) 07900 246529

**North West / North Wales / West Midlands**  
Andy Hargreaves  
[ahargreaves@nsk-uk.com](mailto:ahargreaves@nsk-uk.com) 07872 824175

**Yorkshire & North East**  
Jo Dronfield  
[jdrnfield@nsk-uk.com](mailto:jdrnfield@nsk-uk.com) 07545 786619

**East Midlands**  
Graham Brown  
[gbrown@nsk-uk.com](mailto:gbrown@nsk-uk.com) 07541 864641

**South Central / London SW, W, NW**  
Steve Barnes  
[sbarnes@nsk-uk.com](mailto:sbarnes@nsk-uk.com) 07889 181099

**South East & East Anglia**  
Francois Faro  
[ffaro@nsk-uk.com](mailto:ffaro@nsk-uk.com) 07900 246822

**Field Service Engineer (North)**  
Jamie Duncombe  
[jduncombe@nsk-uk.com](mailto:jduncombe@nsk-uk.com) 07780 105457

**Field Service Engineer (South)**  
Lee Morris  
[lmorris@nsk-uk.com](mailto:lmorris@nsk-uk.com) 07435 783707

**NSK United Kingdom LTD** [www.mynsk.co.uk](http://www.mynsk.co.uk)

Office 4, Gateway 1000 Arlington Business Park, Whittle Way Stevenage, SG1 2FP, UK - [info@nsk-uk.com](mailto:info@nsk-uk.com) - 0800 6341909



\*Offers cannot be combined. NSK reserve the right to change prices and offers without notice. Published prices are RRP's and exclude VAT and installation. Offers only available through participating dealers. Additional warranties available from NSK UK on new products and NSK UK repairs. To qualify for specific offers please send a copy/invoice to NSK UK at Office 4, Gateway 1000, Arlington Business Park, Whittle Way, Stevenage, Hertfordshire SG1 2FP or email [sales@nsk-uk.com](mailto:sales@nsk-uk.com). Offers end 31/08/23.